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```
# axis marking for wavelength
wavelength_start = spectrometer_offset
wavelength_end = spectrometer_offset + 3648*spectrometer_resolution

plot.delete("wavelength")

for x in range(100,1200,50):
    if (x > wavelength_start) and (x < wavelength_end):
        # print(x)
        # print(wavelength_to_pixel(x))

        plot.create_line(PLOT_OFFSETX +
            wavelength_to_pixel(x)*plot_width/PIXEL_NUMBER, plot_height +
            PLOT_OFFSETY, PLOT_OFFSETX + \
                wavelength_to_pixel(x)*plot_width/PIXEL_NUMBER,
                plot_height + PLOT_OFFSETY + 10,
                tags="wavelength")
        plot.create_text(PLOT_OFFSETX +
            wavelength_to_pixel(x)*plot_width/PIXEL_NUMBER, plot_height +
            PLOT_OFFSETY + 20, font=("Arial Bold",
            18),text=str(x),fill='black', \
                tags="wavelength")

spectrometer_offset_output = "Spectrometer offset: " +
    str(spectrometer_offset) + str(" nm")
spectrometer_resolution_output = "Spectrometer resolution: " +
    str(spectrometer_resolution) + str(" nm/pixel")
wavelength_output = "Wavelength in nm          " + spectrometer_offset_output
    + ",          " + spectrometer_resolution_output + ")"
plot.create_text(PLOT_OFFSETX, plot_height + PLOT_OFFSETY + 45,
    font=("Arial Bold",18),text=wavelength_output,fill='black',
    tags="wavelength", anchor="w")

# calibrate the spectrometer using the two choosen peaks from the neon
# spectrum
def calibrate_sensor():

    global spectrometer_offset, spectrometer_resolution

    # get the wavelengths of the two choosen peaks
    wavelength_1 = peaks[int(peak_num_entry.get())-1]
    wavelength_2 = peaks[int(peak_num_entry2.get())-1]

    # calculate the spectrometer resolution in nm/pixel
    spectrometer_resolution = round(((703.24 - 585.25) / (wavelength_2 -
        wavelength_1),4)
    # calculate the offset of the spectrometer, which is defined as the
    # wavelength at the first pixel
    spectrometer_offset = round(585.25 - wavelength_1 *
        spectrometer_resolution,2)

    output_text.insert(END, "\nCalibration successful!\n\n")
    output_text.insert(END, "Spectrometer resolution: " +
        str(spectrometer_resolution) + "nm/pixel")
    output_text.insert(END, "\nSpectrometer offset: " +
```

```
    str(spectrometer_offset) + "nm")

is_calibrated = True
Plot_Sensordata_Wavelength()

# writing configuration file
file = open("EasyCalibrate_config.txt", "w")
file.write("spectrometer_offset " + str(spectrometer_offset) + "\n")
file.write("spectrometer_resolution " + str(spectrometer_resolution))
file.close()

def wavelength_to_pixel(wavelength):

    global spectrometer_offset, spectrometer_resolution

    # print(wavelength)
    pixel = (wavelength - spectrometer_offset) / spectrometer_resolution

    # print(pixel)
    return pixel

def close_peak_window():

    plot.delete("peaks")
    peak_window.destroy()

def find_and_show_peaks():

    global peaks, peak_num_entry, peak_num_entry2, output_text, peak_window

    # Receive sensor values
    sensor_data = [pointer[pixel] for pixel in range(PIXEL_NUMBER)]

    # Finding peaks with the find_peaks function of SciPy
    peaks, _ = find_peaks(sensor_data, height=10000, distance=50)

    # Create a new window for the output of the peaks found
    peak_window = Toplevel(master)
    peak_window.title("Found peaks")
    peak_window.geometry("600x600")

    # Text field for displaying the peaks found
    output_text = Text(peak_window, font=("Arial", 12), width=60, height=20)
    output_text.grid(row=0, column=0, columnspan=2)

    # Entry fields for peak numbers
    peak_num_label = Label(peak_window, text="Which peak number for
        585,25nm?", font=("Arial", 12))
    peak_num_label.grid(row=1, column=0)
    peak_num_entry = Entry(peak_window, font=("Arial", 12), width=5)
    peak_num_entry.grid(row=1, column=1)

    peak_num_label2 = Label(peak_window, text="Which peak number for
        703,24nm?", font=("Arial", 12))
    peak_num_label2.grid(row=2, column=0)
    peak_num_entry2 = Entry(peak_window, font=("Arial", 12), width=5)
```

```
peak_num_entry2.grid(row=2, column=1)

# Calibration button
calibrate = Button(peak_window, text="Calibrate", font=("Arial", 12),
    command=calibrate_sensor)
calibrate.grid(row=3, column=0, columnspan=2)

# Button for closing the window
exit_button = Button(peak_window, text="Exit", font=("Arial", 12),
    command=close_peak_window)
exit_button.grid(row=4, column=0, columnspan=2)

# Output of the peaks found in the text field
output_text.insert(END, "Found peaks:\n\nNumber      Pixel#      Counts\nWavelength [nm]\n")
for idx, peak_index in enumerate(peaks):
    wavelength = spectrometer_offset + peak_index *
        spectrometer_resolution
    if wavelength > 0:
        output_text.insert(END, f"{idx+1:>8}:          {peak_index:>5}\n
            {sensor_data[peak_index]:>5}          {wavelength:6.2f}\n")
    else:
        output_text.insert(END, f"{idx+1:>8}:          {peak_index:>5}\n
            {sensor_data[peak_index]:>5}          -\n")

# Update plot and number peaks consecutively
for idx, peak_index in enumerate(peaks):
    counts = sensor_data[peak_index]
    x = int(peak_index / PIXEL_NUMBER * plot_width)
    y = int(counts / COUNT_MAX * plot_height)

    # Update plotting of the peak marker above the PeakPlot and number
    # peaks consecutively
    plot.create_text(x + PLOT_OFFSETX, plot_height + PLOT_OFFSETY - y -
        10,
                    font=("Arial", 10), text=f"Peak {idx+1}",
                    fill="red", tags="peaks")

print("EasyCalibration V1.1\nSearching for camera: ")

# open external DLL
libe9u = ctypes.WinDLL('./libe9u_LSMD_x64.dll')

# define argument and return types for the used functions
libe9u.e9u_LSMD_search_for_camera.argtype = ctypes.c_uint
libe9u.e9u_LSMD_search_for_camera.restype = ctypes.c_int

libe9u.e9u_LSMD_start_camera_async.argtype = ctypes.c_uint
libe9u.e9u_LSMD_start_camera_async.restype = ctypes.c_int

libe9u.e9u_LSMD_set_times_us.argtypes = (ctypes.c_uint, ctypes.c_uint,
    ctypes.c_uint)
libe9u.e9u_LSMD_set_times_us.restype = ctypes.c_int

libe9u.e9u_LSMD_get_next_frame.argtype = ctypes.c_uint
libe9u.e9u_LSMD_get_next_frame.restype = ctypes.c_int
```

```
libe9u.e9u_LSMD_get_pixel_pointer.argtypes = (ctypes.c_uint, ctypes.c_uint)
libe9u.e9u_LSMD_get_pixel_pointer.restype = ctypes.POINTER(ctypes.c_uint16)

# Search for a suitable camera on all USB ports and quit with returning the
  error code, if no camera is found
i_status = libe9u.e9u_LSMD_search_for_camera(0)
if i_status != 0:
    print("No camera found! Error Code: " + str(i_status))
    exit(1)

print("Starting camera: ", end='')
libe9u.e9u_LSMD_start_camera_async(0)

# getting the pointer to the array containing the sensor data
pointer = libe9u.e9u_LSMD_get_pixel_pointer(0, 0)

# defining the master window for graphical output
master = Tk()

# getting the size of the master window
screen_width = master.winfo_screenwidth()
screen_height = master.winfo_screenheight()

# setting the size of the display window to cover nearly the complete screen
master.geometry(str(screen_width - 50) + "x" + str(screen_height - 100) +
    "+10+20")

# defining the dimensions for the plot area
plot_width = screen_width - 150
plot_height = screen_height - 300

# defining and packing the control/output widgets
statusline = Frame(master)
statusline.pack(side='top')
plot = Canvas(master)
plot.pack(side='bottom', fill=BOTH, expand=YES)

# output label for program name and version number
output_peak_width = Label(statusline, text="EasyCalibrate V1.1 /
    www.eureca.de\nx-axis: pixel number/wavelength\ny-axis: counts",
    justify=LEFT, font=("Arial Bold", 18))
output_peak_width.pack(side='left', padx=0)

# defining slider for integration time and setting it to 10
faktor = IntVar()
slider_exp_time = Scale(statusline, from_=1, to=1000, length=plot_width/3,
    orient=HORIZONTAL, label="Integration time:", font=("Arial Bold", 18))
slider_exp_time.pack(side='left', padx=50)
slider_exp_time.set(100)

# defining two radio buttons for switching the integration time between µs
  and ms
Radiobutton_us = Radiobutton(statusline, text="µs", font=("Arial Bold",
    18), variable=faktor, value=1)
Radiobutton_us.pack(side='left', padx=20)
```

```
Radiobutton_us.invoke()
Radiobutton_ms = Radiobutton(statusline, text="ms", font=("Arial Bold",
    18), variable=faktor, value=1000)
Radiobutton_ms.pack(side='left', padx=20)

start_stop_button = Button(statusline, text="Start/Stop", font=("Arial
    Bold", 18), command=toggle_reading)
start_stop_button.pack(side='left', padx=50)

find_peaks_button = Button(statusline, text="Find Peaks", font=("Arial
    Bold", 18), command=find_and_show_peaks)
find_peaks_button.pack(side='left', padx=20)

# defining the exit button with the closing function
def close_window():
    master.destroy()
exit_button = Button(statusline, text="Exit", font=("Arial Bold", 18),
    command=close_window)
exit_button.pack(side='left', padx=20)

# drawing a rectangular frame for the sensor data
Sensor_Plot = plot.create_rectangle(PLOT_OFFSETX, PLOT_OFFSETY, plot_width
    + PLOT_OFFSETX, plot_height + PLOT_OFFSETY, fill="#fffbe6")

# x-axis marking with pixel number
for x in range(15):
    plot.create_line(PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY, PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY - 10)
    plot.create_text(PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY - 20,font=("Arial Bold", 18),text=int(x*250),fill='black')

    plot.create_line(PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY, PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY + plot_height, fill="#dddddd")

# y-axis marking with count number
for y in range(14):
    plot.create_line(PLOT_OFFSETX-10, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX,
        PLOT_OFFSETX, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX)
    text_output = str(y*5) + "k"
    plot.create_text(PLOT_OFFSETX-40, plot_height + PLOT_OFFSETY -
        y*5000*plot_height/COUNT_MAX,font=("Arial Bold",
        18),text=text_output,fill='black')

    plot.create_line(PLOT_OFFSETX + plot_width, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX,
        PLOT_OFFSETX, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX, fill="#dddddd")

spectrometer_offset = 0.0
spectrometer_resolution = 0.0

# check if there is a cofiguration file in the same directory
```

```
if os.path.exists("EasyCalibrate_config.txt"):
    # reading configuration file, values indicated here will overwrite the
    # default values
    file = open("EasyCalibrate_config.txt", "r")
    for line in file:
        config_variable = line.split()[0]
        config_value = line.split()[1]
        # print(str(config_variable) + " = " + str(config_value))
        # if config_variable == "spectrometer_name":
        #     spectrometer_name = config_value
        if config_variable == "spectrometer_offset":
            spectrometer_offset = float(config_value)
        if config_variable == "spectrometer_resolution":
            spectrometer_resolution = float(config_value)

    file.close()

    # print the sensor wavelength at the lower side of the plot window
    Plot_Sensordata_Wavelength()
else:
    plot.create_text(PLOT_OFFSETX + 100, plot_height + PLOT_OFFSETY + 20,
                    font=("Arial Bold", 18),text="not calibrated",fill='black',
                    tags="wavelength")

def update_plot():
    if not master.wininfo_exists():
        return

    if is_reading:
        # getting the integration time from the slider and using this value for
        # the exposure time as well as for the frame time
        exposure_time = int(slider_exp_time.get()) * faktor.get()
        frame_time = exposure_time;

        # reading out the camera; for details refer to the EasyAccess
        # documentation
        libe9u.e9u_LSMD_set_times_us (0, exposure_time, frame_time)
        libe9u.e9u_LSMD_get_next_frame (0)

        x_old = 0
        y_old = 0

        # deleting the old data points
        plot.delete("data")

        for pixel in range(0, PIXEL_NUMBER):

            counts = pointer[pixel]

            # scaling the sensor data to fit into the plot region
            x = int(pixel / PIXEL_NUMBER * plot_width)
            y = int(counts / COUNT_MAX * plot_height)

            # plotting the data point via connecting the current data with the
            # last one
            plot.create_line(x_old + PLOT_OFFSETX,plot_height + PLOT_OFFSETY -
```

```
        y_old,x + PLOT_OFFSETX,plot_height + PLOT_OFFSETY - y,  
        tags="data")  
  
        x_old = x  
        y_old = y  
  
    plot.update()  
  
    master.after(1, update_plot)  
  
update_plot()  
  
master.mainloop()
```

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