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4 EasyCalibrate Python-Code

```
# EasyCalibrate.py V1.1
#
# Python display tool for line scan cameras by EURECA Messtechnik GmbH
# - detects the camera via the EasyAccess DLL
# - reads out the camera and displays the recorded sensor data
# - integration time can be adjusted between 1µs and 1000ms
#
# For details please refer to: www.eureca.de

# import library for widgets
from tkinter import *

import os
import numpy as np
from scipy.signal import find_peaks

# import library for handling external camera DLL
import ctypes

# basic values of the used linear sensor; e.g. for e9u-LSMD-TCD1304-STD:
    3648 pixel; 16bit data
PIXEL_NUMBER = int(3648)
COUNT_MAX = int(65536)

# offsets for the plot window; needed for displaying additional data
PLOT_OFFSETX = 70
PLOT_OFFSETY = 50

# variable to test the camera status and the wavelength calibration
is_reading = True

# Function for starting/stopping the sensor readout
def toggle_reading():
    global is_reading
    is_reading = not is_reading

def Plot_Sensordata_Wavelength():

    # axis marking for wavelength
    wavelength_start = spectrometer_offset
    wavelength_end = spectrometer_offset + 3648*spectrometer_resolution

    plot.delete("wavelength")

    for x in range(100,1200,50):
        if (x > wavelength_start) and (x < wavelength_end):
            # print(x)
            # print(wavelength_to_pixel(x))

            plot.create_line(PLOT_OFFSETX +
                wavelength_to_pixel(x)*plot_width/PIXEL_NUMBER, plot_height +
                PLOT_OFFSETY, PLOT_OFFSETX + \
                    wavelength_to_pixel(x)*plot_width/PIXEL_NUMBER,
                plot_height + PLOT_OFFSETY + 10,
```

```
        tags="wavelength")
    plot.create_text(PLOT_OFFSETX +
        wavelength_to_pixel(x)*plot_width/PIXEL_NUMBER, plot_height +
        PLOT_OFFSETY + 20, font=("Arial Bold",
        18),text=str(x),fill='black', \
        tags="wavelength")

spectrometer_offset_output = "Spectrometer offset: " +
    str(spectrometer_offset) + str(" nm")
spectrometer_resolution_output = "Spectrometer resolution: " +
    str(spectrometer_resolution) + str(" nm/pixel")
wavelength_output = "Wavelength in nm          " + spectrometer_offset_output
    + ",          " + spectrometer_resolution_output + ")"
plot.create_text(PLOT_OFFSETX, plot_height + PLOT_OFFSETY + 45,
    font=("Arial Bold",18),text=wavelength_output,fill='black',
    tags="wavelength", anchor="w")

# calibrate the spectrometer using the two choosen peaks from the neon
# spectrum
def calibrate_sensor():

    global spectrometer_offset, spectrometer_resolution

    # get the wavelengths of the two choosen peaks
    wavelength_1 = peaks[int(peak_num_entry.get())-1]
    wavelength_2 = peaks[int(peak_num_entry2.get())-1]

    # calculate the spectrometer resolution in nm/pixel
    spectrometer_resolution = round((703.24 - 585.25) / (wavelength_2 -
        wavelength_1),4)
    # calculate the offset of the spectrometer, which is defined as the
    # wavelength at the first pixel
    spectrometer_offset = round(585.25 - wavelength_1 *
        spectrometer_resolution,2)

    output_text.insert(END, "\nCalibration successful!\n\n")
    output_text.insert(END, "Spectrometer resolution: " +
        str(spectrometer_resolution) + "nm/pixel")
    output_text.insert(END, "\nSpectrometer offset: " +
        str(spectrometer_offset)+"nm")

    is_calibrated = True
    Plot_Sensordata_Wavelength()

    # writing configuration file
    file = open("EasyCalibrate_config.txt", "w")
    file.write("spectrometer_offset " + str(spectrometer_offset) + "\n")
    file.write("spectrometer_resolution " + str(spectrometer_resolution))
    file.close()

def wavelength_to_pixel(wavelength):

    global spectrometer_offset, spectrometer_resolution

    # print(wavelength)
    pixel = (wavelength - spectrometer_offset) / spectrometer_resolution
```

```
# print(pixel)
return pixel

def close_peak_window():

    plot.delete("peaks")
    peak_window.destroy()

def find_and_show_peaks():

    global peaks, peak_num_entry, peak_num_entry2, output_text, peak_window

    # Receive sensor values
    sensor_data = [pointer[pixel] for pixel in range(PIXEL_NUMBER)]

    # Finding peaks with the find_peaks function of SciPy
    peaks, _ = find_peaks(sensor_data, height=10000, distance=50)

    # Create a new window for the output of the peaks found
    peak_window = Toplevel(master)
    peak_window.title("Found peaks")
    peak_window.geometry("600x600")

    # Text field for displaying the peaks found
    output_text = Text(peak_window, font=("Arial", 12), width=60, height=20)
    output_text.grid(row=0, column=0, columnspan=2)

    # Entry fields for peak numbers
    peak_num_label = Label(peak_window, text="Which peak number for
        585,25nm?", font=("Arial", 12))
    peak_num_label.grid(row=1, column=0)
    peak_num_entry = Entry(peak_window, font=("Arial", 12), width=5)
    peak_num_entry.grid(row=1, column=1)

    peak_num_label2 = Label(peak_window, text="Which peak number for
        703,24nm?", font=("Arial", 12))
    peak_num_label2.grid(row=2, column=0)
    peak_num_entry2 = Entry(peak_window, font=("Arial", 12), width=5)
    peak_num_entry2.grid(row=2, column=1)

    # Calibration button
    calibrate = Button(peak_window, text="Calibrate", font=("Arial", 12),
        command=calibrate_sensor)
    calibrate.grid(row=3, column=0, columnspan=2)

    # Button for closing the window
    exit_button = Button(peak_window, text="Exit", font=("Arial", 12),
        command=close_peak_window)
    exit_button.grid(row=4, column=0, columnspan=2)

    # Output of the peaks found in the text field
    output_text.insert(END, "Found peaks:\n\nNumber      Pixel#      Counts\nWavelength [nm]\n")
    for idx, peak_index in enumerate(peaks):
        wavelength = spectrometer_offset + peak_index *
```

```
spectrometer_resolution
if wavelength > 0:
    output_text.insert(END, f"{idx+1:>8}:           {peak_index:>5}
                        {sensor_data[peak_index]:>5}       {wavelength:6.2f}\n")
else:
    output_text.insert(END, f"{idx+1:>8}:           {peak_index:>5}
                        {sensor_data[peak_index]:>5}       -\n")

# Update plot and number peaks consecutively
for idx, peak_index in enumerate(peaks):
    counts = sensor_data[peak_index]
    x = int(peak_index / PIXEL_NUMBER * plot_width)
    y = int(counts / COUNT_MAX * plot_height)

# Update plotting of the peak marker above the PeakPlot and number
# peaks consecutively
plot.create_text(x + PLOT_OFFSETX, plot_height + PLOT_OFFSETY - y -
                10,
                font=("Arial", 10), text=f"Peak {idx+1}",
                fill="red", tags="peaks")

print("EasyCalibration V1.1\nSearching for camera: ")

# open external DLL
libe9u = ctypes.WinDLL('./libe9u_LSMD_x64.dll')

# define argument and return types for the used functions
libe9u.e9u_LSMD_search_for_camera.argtype = ctypes.c_uint
libe9u.e9u_LSMD_search_for_camera.restype = ctypes.c_int

libe9u.e9u_LSMD_start_camera_async.argtype = ctypes.c_uint
libe9u.e9u_LSMD_start_camera_async.restype = ctypes.c_int

libe9u.e9u_LSMD_set_times_us.argtypes = (ctypes.c_uint, ctypes.c_uint,
                                         ctypes.c_uint)
libe9u.e9u_LSMD_set_times_us.restype = ctypes.c_int

libe9u.e9u_LSMD_get_next_frame.argtype = ctypes.c_uint
libe9u.e9u_LSMD_get_next_frame.restype = ctypes.c_int

libe9u.e9u_LSMD_get_pixel_pointer.argtypes = (ctypes.c_uint, ctypes.c_uint)
libe9u.e9u_LSMD_get_pixel_pointer.restype = ctypes.POINTER(ctypes.c_uint16)

# Search for a suitable camera on all USB ports and quit with returning the
# error code, if no camera is found
i_status = libe9u.e9u_LSMD_search_for_camera(0)
if i_status != 0:
    print("No camera found! Error Code: " + str(i_status))
    exit(1)

print("Starting camera: ", end='')
libe9u.e9u_LSMD_start_camera_async(0)

# getting the pointer to the array containing the sensor data
pointer = libe9u.e9u_LSMD_get_pixel_pointer(0, 0)
```

```
# defining the master window for graphical output
master = Tk()

# getting the size of the master window
screen_width = master.winfo_screenwidth()
screen_height = master.winfo_screenheight()

# setting the size of the display window to cover nearly the complete screen
master.geometry(str(screen_width - 50) + "x" + str(screen_height - 100) +
    "+10+20")

# defining the dimensions for the plot area
plot_width = screen_width - 150
plot_height = screen_height - 300

# defining and packing the control/output widgets
statusline = Frame(master)
statusline.pack(side='top')
plot = Canvas(master)
plot.pack(side='bottom', fill=BOTH, expand=YES)

# output label for program name and version number
output_peak_width = Label(statusline, text="EasyCalibrate V1.1 /
    www.eureca.de\nx-axis: pixel number/wavelength\ny-axis: counts",
    justify=LEFT, font=("Arial Bold", 18))
output_peak_width.pack(side='left', padx=0)

# defining slider for integration time and setting it to 10
faktor = IntVar()
slider_exp_time = Scale(statusline, from_=1, to=1000, length=plot_width/3,
    orient=HORIZONTAL, label="Integration time:", font=("Arial Bold", 18))
slider_exp_time.pack(side='left', padx=50)
slider_exp_time.set(100)

# defining two radio buttons for switching the integration time between µs
and ms
Radiobutton_us = Radiobutton(statusline, text="µs", font=("Arial Bold",
    18), variable=faktor, value=1)
Radiobutton_us.pack(side='left', padx=20)
Radiobutton_us.invoke()
Radiobutton_ms = Radiobutton(statusline, text="ms", font=("Arial Bold",
    18), variable=faktor, value=1000)
Radiobutton_ms.pack(side='left', padx=20)

start_stop_button = Button(statusline, text="Start/Stop", font=("Arial
    Bold", 18), command=toggle_reading)
start_stop_button.pack(side='left', padx=50)

find_peaks_button = Button(statusline, text="Find Peaks", font=("Arial
    Bold", 18), command=find_and_show_peaks)
find_peaks_button.pack(side='left', padx=20)

# defining the exit button with the closing function
def close_window():
    master.destroy()
exit_button = Button(statusline, text="Exit", font=("Arial Bold", 18),
```

```
command=close_window)
exit_button.pack(side='left', padx=20)

# drawing a rectangular frame for the sensor data
Sensor_Plot = plot.create_rectangle(PLOT_OFFSETX, PLOT_OFFSETY, plot_width
    + PLOT_OFFSETX, plot_height + PLOT_OFFSETY, fill="#fffbe6")

# x-axis marking with pixel number
for x in range(15):
    plot.create_line(PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY, PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY - 10)
    plot.create_text(PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY - 20,font=("Arial Bold", 18),text=int(x*250),fill='black')

    plot.create_line(PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY, PLOT_OFFSETX + x*250*plot_width/PIXEL_NUMBER,
        PLOT_OFFSETY + plot_height, fill="#dddddd")

# y-axis marking with count number
for y in range(14):
    plot.create_line(PLOT_OFFSETX-10, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX,
        PLOT_OFFSETX, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX)
    text_output = str(y*5) + "k"
    plot.create_text(PLOT_OFFSETX-40, plot_height + PLOT_OFFSETY -
        y*5000*plot_height/COUNT_MAX,font=("Arial Bold",
        18),text=text_output,fill='black')

    plot.create_line(PLOT_OFFSETX + plot_width, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX,
        PLOT_OFFSETX, PLOT_OFFSETY + plot_height -
        y*5000*plot_height/COUNT_MAX, fill="#dddddd")

spectrometer_offset = 0.0
spectrometer_resolution = 0.0

# check if there is a configuration file in the same directory
if os.path.exists("EasyCalibrate_config.txt"):
    # reading configuration file, values indicated here will overwrite the
    # default values
    file = open("EasyCalibrate_config.txt", "r")
    for line in file:
        config_variable = line.split()[0]
        config_value = line.split()[1]
        # print(str(config_variable) + " = " + str(config_value))
        # if config_variable == "spectrometer_name":
        #     spectrometer_name = config_value
        if config_variable == "spectrometer_offset":
            spectrometer_offset = float(config_value)
        if config_variable == "spectrometer_resolution":
            spectrometer_resolution = float(config_value)

    file.close()
```


```
# print the sensor wavelength at the lower side of the plot window
Plot_Sensordata_Wavelength()
else:
    plot.create_text(PLOT_OFFSETX + 100, plot_height + PLOT_OFFSETY + 20,
        font=("Arial Bold", 18),text="not calibrated",fill='black',
        tags="wavelength")

def update_plot():
    if not master.wininfo_exists():
        return

    if is_reading:
        # getting the integration time from the slider and using this value for
        # the exposure time as well as for the frame time
        exposure_time = int(slider_exp_time.get()) * faktor.get()
        frame_time = exposure_time;

        # reading out the camera; for details refer to the EasyAccess
        # documentation
        libe9u.e9u_LSMD_set_times_us (0, exposure_time, frame_time)
        libe9u.e9u_LSMD_get_next_frame (0)

        x_old = 0
        y_old = 0

        # deleting the old data points
        plot.delete("data")

        for pixel in range(0, PIXEL_NUMBER):

            counts = pointer[pixel]

            # scaling the sensor data to fit into the plot region
            x = int(pixel / PIXEL_NUMBER * plot_width)
            y = int(counts / COUNT_MAX * plot_height)

            # plotting the data point via connecting the current data with the
            # last one
            plot.create_line(x_old + PLOT_OFFSETX,plot_height + PLOT_OFFSETY -
                y_old,x + PLOT_OFFSETX,plot_height + PLOT_OFFSETY - y,
                tags="data")

            x_old = x
            y_old = y

        plot.update()

        master.after(1, update_plot)

update_plot()

master.mainloop()
```

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